







HOW TO USE BD™ ESwab™ and BD™ Universal Viral Transport



BD Life Sciences
Diagnostic Systems
7 Loveton Circle
Sparks, MD 21152
Tel: 800.638.8663
www.bd.com/ds

UNC Microbiology Lab

101 Manning Drive • Chapel Hill, NC 27514
984.974.1805

<http://www.unccmedicalcenter.org/mclendon-clinical-laboratories/>



Lawson # 309762

BD Cat No.:220245

BD ESwab Regular Collection Kit

White cap w/liquid amies and a regular flocced swab

For the collection and transport of clinical specimens containing aerobes, anaerobes and fastidious bacteria

Anaerobic | Aerobic wounds, Eye Cultures, CF Sputum (deep pharyngeal), MRSA, VRE, yeast screens

Lawson # 143675

BD Cat No.:220532

BD ESwab Flexible Minitip Kit

Blue cap w/liquid amies & a flex minitip flocced swab

For the collection and transport of clinical specimens containing aerobes, anaerobes and fastidious bacteria

Ear culture, NP bacterial culture, Sinus culture



Lawson # 050750

BD Cat No.:220531

BD™ Universal Viral Transport

with flexible mini-tip flocced swab

For the collection and transport of clinical specimens containing viruses.

Flu, Flu | RSV, RVP, VZV, Adenovirus PCR (Eye, Peri-Rectal Swab)

HSV PCR (Lesion Swabs)

VZV PCR (Lesion Swabs)

Ureaplasma PCR







Mycoplasma PCR

Chlamydia trach. culture (NP or throat swab)

Chlamydia pneumonia PCR (NP or throat swab)



HOW TO USE BD ESwab And BD Universal Viral Transport



BD Life Sciences
Diagnostic Systems
7 Loveton Circle
Sparks, MD 21152
Tel: 800.638.8663



Lawson # 050440

Cat. No. 413

HealthLink Liquid Stuart Swab

Orange Cap with regular aluminum wire

For the collection and transport of clinical specimens containing aerobes.

Bordetella pertussis

Bordetella parapertussis



Lawson # 050439

Cat No.: 4019

HealthLink Liquid Stuart DoubleSwab

Red Cap with two (2) swabs

For the collection and transport of clinical specimens containing aerobes.

Rapid Strep

Group A Strep Culture

Group B Strep Culture

Wet Prep, Genital

Fungal Culture

UNC Microbiology Lab

101 Manning Drive • Chapel Hill, NC 27514

984.974.1805

<http://www.uncm edicalcenter.org/m clendon-clinical-laboratories/>